



ELSEVIER

Analytica Chimica Acta 403 (2000) 249–258

ANALYTICA
CHIMICA
ACTA

www.elsevier.com/locate/aca

Synthetic zeolites as sorbent material for on-line preconcentration of copper traces and its determination using flame atomic absorption spectrometry

Yaneira Petit de Peña^{a,*}, Wilmer López^a, José Luis Burguera^a, Marcela Burguera^a,
Máximo Gallignani^a, Rosario Brunetto^a, Pablo Carrero^a, Carlos Rondon^a, Freddy Imbert^b

^a IVAIQUIM (Venezuelan Andean Institute for Chemical Research), Faculty of Sciences, La Hechicera, University of Los Andes, Mérida 5101-A, Venezuela

^b Laboratory of Kinetics and Catalysis, Faculty of Sciences, University of Los Andes, Mérida 5101-A, Venezuela

Received 24 March 1999; received in revised form 20 July 1999; accepted 26 July 1999

Abstract

This work assesses for the first time the potential of synthetic zeolites as adsorptive material for the preconcentration of copper traces. This involves the formation of a copper chelate, its adsorption on sodium lynde Type A (Na-LTA) and sodium faujasite Type X (Na-Fau) synthetic zeolite minicolumns, and its elution with 300 μ l of methyl isobutyl ketone. Using 5–80 mg of either sorbents materials with a loading time of 2 min per sample preconcentration factors of copper ranging from 35–125 for Na-LTA and 30–65 for Na-Fau were readily achieved related to the direct introduction of aqueous solutions into an atomic absorption spectrometer. To evaluate the performance of these materials, a flow injection system was used for samples of water from different sources, and the accuracy was assessed through recovery experiments and the analysis of certified reference materials. A parallel experiment was performed with polygosyl bonded silica reversed-phase sorbent with octadecyl functional groups (RP-C₁₈) as sorbent in the minicolumn in order to study the performance of this well-known material as a reference for the new sorbents. The 3σ detection limits, relative standard deviation and linear calibration graphs were, respectively, 0.1 ng ml⁻¹, 2.6% and 2.5–30 ng ml⁻¹ for Na-LTA and 0.4 ng ml⁻¹, 2.8% and 2.5–40 ng ml⁻¹ for Na-Fau for preconcentration times of 2 min. ©2000 Elsevier Science B.V. All rights reserved.

Keywords: On-line preconcentration; Synthetic zeolites; Flow injection; Copper; Atomic absorption spectroscopy

1. Introduction

The papers on flow injection atomic spectrometry (FI-AS) systems used for preconcentration purposes that have been found in the literature could be classi-

fied in three groups: liquid–liquid extraction, column preconcentration and precipitation, according to the basic mechanism of the preconcentration step [1,2].

Particularly, solid-phase extraction (SPE) is a preconcentration technique of rapidly growing importance in trace metal determination, which can be easily adapted to FI on line separation and preconcentration systems [3]. The main advantage of this is the possibility of using a relatively simple detection

* Corresponding author. Tel.: +58-74-401375; fax: +58-74-401286

E-mail address: ypetit@ciens.ula.ve (Y.P. de Peña)